

E-PF-P

Advanced Pipette Filler

Features

- Digital pipette filler compatible with glass or plastic serological pipettes (0.1 -100mL)
- Intuitive LCD screen displays battery status and pipetting speed
- Convenient and precise, 6 speed levels for aspiration and dispensing
- Ergonomic lightweight design for fatigue-free and efficient pipetting
- The quiet yet powerful motor with pump fills a 25mL pipette in less than 5 seconds on level 6
- Lithium-ion battery
- Includes a charger and stand




E-PF-M

Manual Pipette Filler

Features

- Economic and easy to use
- Compatible with glass or plastic serological pipettes (0.1 - 100mL)
- 3µm replaceable hydrophobic filter
- Easy to clean and maintain
- Aspiration and dispensing comfortably achieved with a sensitive lever
- Convenient button to expel residual liquid when using blow-out pipette
- Manually compress the suction bellow before attaching the pipette

Autoclavable Silicone Adapters for 0.1-100ml pipettes - replace every 6 months



Accessories

Code	E-PF-F3.0	E-PF-MSA	E-PF-F0.45	E-PF-PSA
Picture				
Type	3.0µm Hydrophobic Filter	Silicone Adapter	0.45µm Hydrophobic Filter	Silicone Adapter
Included	1pc with E-PF-M	1pcs with E-PF-M	2pcs with E-PF-P	1pcs with E-PF-P

Specifications		E-PF-M	E-PF-P
Aspiration speed		Manual control, variable level	25mL <5s, 6 Levels
Dispensing method		Gravity	Gravity / Motor, 6 Levels
Pipette Types		Glass or Plastic Pipette (0.1-100mL), Pasteur Pipettes	Glass or Plastic Pipette (0.1-100mL), Pasteur Pipettes
Battery	Type	-	Lithium-ion
	Full charge	-	2-3hrs
	Run time	-	>8hrs of intermittent use
Hydrophobic filter		E-PF-F3.0 (3.0µm)	E-PF-F0.45 (0.45µm)
Silicone adapter		E-PF-MSA	E-PF-PSA
Weight		120g	200g